



THIS DATA SHEET IS FOR IDENTIFICATION PURPOSES ONLY. NOMINAL DIMENSIONS IN MM.

CABLE TYPE C2042

TECHNICAL DATA:

TERMINATION:-

NOMINAL OPERATING VOLTAGE 210KV DC
- USING AN ESSEX X-RAY APPROVED RECEPTACLE
CABLE LENGTH TOLERANCE +/- 100mm

CABLE:- VOLTAGE WITHSTAND 160KV DC
50KV (RMS) AC

FOR FULL CABLE DATA SEE ESSEX CABLE DATA SHEET

ACCESSORIES:

CLAMPING SLEEVES:-

STANDARD FLANGE T20172/A

QUICK LOCK T20173

INSTALLATION GREASE P00075

INSTALLATION:

WARNING:

CABLE SHOULD ONLY BE FITTED BY A TRAINED PERSON.
OIL MUST NOT BE USED ON THE RUBBER TERMINATION, ONLY SUITABLE HIGH VOLTAGE GREASE.
ENSURE BY ADJUSTING THE CLAMPING SLEEVE THAT WHEN THE TERMINATION IS FULLY INSERTED IN THE SOCKET THERE IS A GAP OF 6mm BETWEEN THE MATING FACES OF THE CLAMPING SLEEVE AND SOCKET BEFORE TIGHTENING THE CLAMPING SLEEVE.

**ESSEX X-RAY &
MEDICAL EQUIPMENT LTD.**

FLITCH INDUSTRIAL ESTATE, CHELMSFORD ROAD
GREAT DUNMOW, ESSEX CM6 1XJ ENGLAND

Telephone 01371 875661

Fax 01371 875665

Email sales@essex-x-ray.com

REVISIONS

NO.	DATE	BY
6	11.05.01	DCN 1455
7	13.07.01	DCN 1534
8	08.07.04	DCN 2047
9	06.03.09	DCN 2905

DUE TO CONTINUING PRODUCT DEVELOPMENT ESSEX X-RAY RESERVES THE RIGHT TO MODIFY DATA WITHOUT NOTICE

DRAWN	APPROVED	DATE	ISSUE	SCALE
JDF		15.09.95	9	N.T.S.

© NOT TO BE REPRODUCED IN WHOLE OR PART WITHOUT PRIOR PERMISSION OF ESSEX X-RAY & MEDICAL EQUIPMENT LTD.

TITLE
210KV SEIFERT RIGHT ANGLE (R30 TYPE)
(CABLE TYPE C2042)

DRAWING No.

DS1006