





TECHNICAL DATA:			ACCESSORIES:				INSTALLATION:					
TERMINATION:— NOMINAL OPERATING VOLTAGE 210KV DC — USING AN ESSEX X-RAY APPROVED RECEPTACLE CABLE LENGTH TOLERANCE +/- 100mm CABLE:— VOLTAGE WITHSTAND 160KV DC 50KV (RMS) AC  FOR FULL CABLE DATA SEE ESSEX CABLE DATA SHEET 			CLAMPING SLEEVES:— STANDARD FLANGE T20172/A QUICK LOCK T20173 INSTALLATION GREASE P00075				WARNING: CABLE SHOULD ONLY BE FITTED BY A TRAINED PERSON. OIL MUST NOT BE USED ON THE RUBBER TERMINATION, ONLY SUITABLE HIGH VOLTAGE GREASE. ENSURE BY ADJUSTING THE CLAMPING SLEEVE THAT WHEN THE TERMINATION IS FULLY INSERTED IN THE SOCKET THERE IS A GAP OF 6mm BETWEEN THE MATING FACES OF THE CLAMPING SLEEVE AND SOCKET BEFORE TIGHTENING THE CLAMPING SLEEVE.					
<div>ESSEX X-RAY & MEDICAL EQUIPMENT LTD.</div> <div>FLITCH INDUSTRIAL ESTATE, CHELMSFORD ROAD GREAT DUNMOW, ESSEX CM6 1XJ ENGLAND Telephone 01371 875661 Fax 01371 875665 Email sales@essex-x-ray.com</div>			REVISIONS			DUE TO CONTINUING PRODUCT DEVELOPMENT ESSEX X-RAY RESERVES THE RIGHT TO MODIFY DATA WITHOUT NOTICE						
			NO.	DATE	BY	DRAWN	APPROVED	DATE	ISSUE	SCALE	© NOT TO BE REPRODUCED IN WHOLE OR PART WITHOUT PRIOR PERMISSION OF ESSEX X-RAY & MEDICAL EQUIPMENT LTD.	
			4	02.01.01	DCN 1405	JDF		15.09.95	7	N.T.S.		
			5	11.05.01	DCN 1455	TITLE 210KV SEIFERT STRAIGHT (R30 TYPE) (CABLE TYPE C2042)						
			6	16.07.01	DCN 1534							
			7	06.03.09	DCN 2905							
					DRAWING No. DS1005							