



High Pressure Dark Matter Detectors

James T White
Texas A&M University
Feb. 26, 2010

DM10 – UCLA

High Pressure Gas?

- Possible Advantages

- Nuclear Recoil Discrimination – how does it compare to liquid phase?
- Room Temperature Operation
 - ease of purification !!!
 - no cryogenics !!!
 - no leveling !!!
- Different targets (Xe, Ar, Ne), ~same approach

- Arguments Against

- Pressure vessel vs cryostat
 - BUT: composites?, pressure held externally?
- Size of vessel
 - density ~ 0.1 g/cc → e.g. $30^{1/3} \sim 3$ x bigger for xenon
- Shielding/surface area
 - Active shield + large water tank – ideal for Room T, HP detector

- Bottom Line

- It ain't over until ...
- AND: future, long term running to collect large samples – may be of use


Outline

- S2 vs S1 discrimination
 - Quick review of previous results
- S1 shape discrimination
 - Recent data
- Plans


S1 vs S2 Discrimination

Xenon – 20 Bar


Neutron Interaction



S2/S1 vs S1



NR @ 10 Bar in 7-PMT cell





$S1$ E_{r-max} 18.5 keV

S2 Projection




30 keV reference

Zoom & Log scale



Note: ~ 5.4 ions/keVNR for $W \sim 22$ eV/ion

Xenon




Neon (+.05%Xenon)

Gammas

Neutrons


$S2/S1$

$S2/S1$



AmBe Source
High energy Neutrons



$S2/S1$




Neutron Beam
 E_r Max ~ 120 KeV

S1 Shape Discrimination


View Region



Windows





100 Bar Test Cell




Argon S1

Zero-field alpha signal?




Neon



AmBe Neutrons
NeXe(0.5%) – S1 vs S2



Pure Neon
Zero Field – S1 Shape








Pure Neon – 68 bar




NeXe(0.5%)

A Fast = area in first 400 ns




NeXe(0.5%)



NeXe(0.5%)

Pulse Shape Discrimination



Conclusions/Plans

- HP noble gases and mixtures have excellent discrimination in both S1 vs S2 and S1 pulse shape discrimination (not so sure about xenon S1 disc. at WIMP scale recoil energies)
- But, with the advanced state of development and higher density of liquid nobles, is there a need for future HP detectors?
- Plans:
 - continue investigation of light yield and properties at low recoil energies using neutron beam
 - look at NR discrimination in xenon up to 0.5 g/cc ($\sim 1.8 \times l_q$)
- With combination: water shield, active shield, thin-walled vessels and appropriate light detection instrumentation, it may turn out that HP detectors have a future for long term WIMP property measurement experiments.

Hard to Choose...



Cryogenic?



Room Temperature?