

ASME BPE-2007
(Revision of ASME BPE-2005)

Bioprocessing Equipment

AN AMERICAN NATIONAL STANDARD



**The American Society of
Mechanical Engineers**

Three Park Avenue • New York, NY 10016

CONTENTS

Foreword	vii
Statements of Policy	viii
Committee Roster	ix
Summary of Changes	xii
Part GR General Requirements	1
GR-1 Introduction	1
GR-2 Scope	1
GR-3 Inspection	1
GR-4 Inspector/Examiner	1
GR-5 Responsibilities	1
GR-6 Access for Inspectors	2
GR-7 Manufacturer's Quality Assurance Program	2
GR-8 Metric	2
GR-9 References	2
GR-10 Terms and Definitions	3
Part SD Design for Sterility and Cleanability	9
SD-1 Introduction	9
SD-2 Scope and Purpose	9
SD-3 General Guidelines	9
SD-4 Specific Guidelines	20
SD-5 Testing and Inspection	54
SD-6 Documentation	55
SD-7 Responsibilities	55
Part DT Dimensions and Tolerances for Stainless Steel Automatic Welding and Hygienic Clamp Tube Fittings and Process Components	56
DT-1 Scope	56
DT-2 Pressure Rating	56
DT-3 Marking	56
DT-4 Materials	56
DT-5 Metal Thickness	56
DT-6 Fitting Dimensions	56
DT-7 Tests	57
DT-8 Tolerances	57
DT-9 Welding Ends	57
DT-10 Hygienic Clamp Ends	57
DT-11 Heat Treatment	57
DT-12 Surface Condition	57
DT-13 Packaging	57
DT-14 Minimum Examination Requirements	57
Part MJ Material Joining	76
MJ-1 Scope	76
MJ-2 Materials	76
MJ-3 Joining Processes and Procedures	76
MJ-4 Weld Joint Design and Preparation	77
MJ-5 Filler Material	77
MJ-6 Weld Acceptance Criteria	78
MJ-7 Inspection, Examination, and Testing	78
MJ-8 Procedure Qualification	84

MJ-9	Performance Qualification	84
MJ-10	Documentation Requirements	84
MJ-11	Passivation	85
Part SF	Stainless Steel and Higher Alloy Product Contact Surface Finishes	86
SF-1	Scope	86
SF-2	Objective	86
SF-3	Applications	86
SF-4	Material	86
SF-5	Inspection and Techniques Employed in the Classification of Product Contact Surface Finishes	86
SF-6	Surface Condition	86
SF-7	Electropolishing Procedure Qualification	87
SF-8	Passivation	87
Part SG	Equipment Seals	89
SG-1	Scope and Purpose	89
SG-2	Seal Classes	89
SG-3	General Provisions for Seals in Bioprocessing Service: User Basic Design Requirement	90
SG-4	Special Provisions for Seals in Bioprocessing Service	99
Part PM	Polymer-Based Materials	103
PM-1	Introduction	103
PM-2	Design Considerations for Polymeric Piping, Tubing, Fittings, Valve Bodies, and Other Components	103
PM-3	Polymer Material Joining	106
PM-4	Polymer Interior Product Contact Surfaces of Piping, Tubing, Fittings, Valve Bodies, and Coated or Lined Vessels	109
PM-5	Materials of Construction	109
Figures		
SD-1	Hygienic Connections	11
SD-2	Nonhygienic Connections	13
SD-3	Flat Gasket Applications	14
SD-4	Recommended and Preferred Drop Designs	18
SD-5	Double Block-and-Bleed Valve Assembly	19
SD-6	Instrument Location Detail: Hygienic Design	20
SD-7	Flexible Hygienic Hose Design	22
SD-8	Tank/Vessel Vent Filters	23
SD-9	Nozzle Design	25
SD-10	Sidewall Instrument Ports	26
SD-11	Dip Tube	27
SD-12	Vessel Design Tangential Nozzles	27
SD-13	Vessel Sight Glass Design	28
SD-14-1	Dip Tube Nozzles: Removable Designs	29
SD-14-2	Side and Bottom Connections	30
SD-15	Agitator Mounting Flanges	31
SD-16	Sight Glass Design	32
SD-17	Internal Support Members	33
SD-18	Mitered Fittings	34
SD-19	Typical Nozzle Detail	34
SD-20	Double Tubesheet Heat Exchanger Bonnet Design	35
SD-21-1	Shaft Coupling Construction	36
SD-21-2	Shaft Coupling Seal Arrangement	37
SD-21-3	Fastener Seal Arrangements	38
SD-21-4	Shaft Steady Bearing	39
SD-21-5	Magnetically-Coupled Mixer (Typical Bottom-Mount)	40
SD-21-6	Double Mechanical Cartridge Seal With Debris Well	42

SD-22-1	Typical Clean Steam System Isometric	43
SD-22-2	Clean Steam Point-of-Use Design	44
SD-22-3	Steam Traps for Clean Steam Systems	45
SD-23	Point-of-Use Piping	46
SD-23-1	Physical Break in Point-of-Use Piping	47
SD-24	Transfer Panel Tolerances	51
SD-25	Transfer Panel Looped Headers	52
SD-26	Transfer Panel Jumpers	53
MJ-1	Acceptable and Unacceptable Weld Profiles for Tube Welds	82
SG-1	Basic Components of a Seal	92
SG-2	Single Dry Running Contacting Seal	93
SG-3	Internally Mounted, Process Lubricated Contact Seal	93
SG-4	Externally Mounted, Process Lubricated Contact Seal	94
SG-5	Double Seal Installation	95
SG-6	Tandem Seal Installation	95
SG-7	Seal Piping and Lubrication Plans	96
SG-8	Gas Lubricated Noncontacting Double Seal	97
SG-9	Tandem Seal With Barrier System	98
SG-10	Typical Packing Installation	98
SG-11	V-Ring Packing for Reciprocating Applications	99
SG-12	Open Cross-Sectional Lip Seal	99
SG-13	Labyrinth Seal	99
SG-14	Typical Angle Valve With Rolling Diaphragm and Orifice	100
SG-15	Example of Sampling Valve With Uniformly Loaded Sliding Seal	101
SG-16	Typical In-Line Diaphragm Valve With Weir	101
SG-17	Typical Ball Valve Configuration	102
PM-1	Acceptable and Unacceptable Weld Profiles for Beadless Welds	108

Tables

SD-1	<i>L/D</i> Dimensions for Flow-Through Tee: Full-Size Standard Straight Tee With Blind Cap	15
SD-2	<i>L/D</i> Dimensions for Flow-Through Tee: Short Outlet Reducing Tee With Blind Cap	16
SD-3	Slope Designations for Gravity-Drained Lines	19
SD-4	Annular Spacing Recommendations for Hygienic Dip Tubes	26
SD-5	Recommended Flow Rates to Achieve 5 fps (1.52 mps)	49
SD-6	Recommended Flow Rates for Cleaning Vertical Cylindrical Vessels Having Dished Heads	49
SD-7	Transfer Panel and Jumper Tolerances	50
DT-1	Nominal O.D. Tubing Sizes	59
DT-2	Hygienic Unions	60
DT-3	Chemical Composition for Automatic Weld Ends, %	60
DT-4	Tangent Lengths	60
DT-5	Final Tolerances for Mechanically Polished Fittings and Process Components	61
DT-5.1	Hygienic Clamp Ferrule Standard Dimensions and Tolerances	62
DT-6	Final Tolerances for Electropolished Fittings and Process Components	63
DT-7	Automatic Tube Weld: 90-deg Elbow	63
DT-8	Automatic Tube Weld: 45-deg Elbow	63
DT-9	Automatic Tube Weld: Straight Tee and Cross	64
DT-10	Automatic Tube Weld: Reducing Tee	64
DT-11	Automatic Tube Weld: Concentric and Eccentric Reducer	65
DT-12	Automatic Tube Weld: Hygienic Clamp Joint, 90-deg Elbow	65
DT-13	Automatic Tube Weld: Hygienic Clamp Joint, 45-deg Elbow	66
DT-14	Automatic Tube Weld: Short Outlet Hygienic Clamp Joint Reducing Tee	66
DT-15	Automatic Tube Weld: Short Outlet Hygienic Clamp Joint Tee	67

DT-16	Hygienic Clamp Joint: 90-deg Elbow	67
DT-17	Hygienic Clamp Joint: 45-deg Elbow	68
DT-18	Hygienic Clamp Joint: Straight Tee and Cross	68
DT-19	Hygienic Clamp Joint: Reducing Tee	69
DT-20	Hygienic Clamp Joint: Short Outlet Reducing Tee	69
DT-21	Hygienic Clamp Joint: Concentric and Eccentric Reducer	70
DT-22	Automatic Tube Weld: Ferrule	71
DT-23	Automatic Tube Weld: 180-deg Return Bend	72
DT-24	Hygienic Clamp Joint: 180-deg Return Bend	72
DT-25	Hygienic Mechanical Joint: Short Outlet Run Tee	73
DT-26	Hygienic Clamp Joint: Tube Weld Concentric and Eccentric Reducer	73
DT-27	Hygienic Clamp Joint: Short Outlet Tee	74
DT-28	Automatic Tube Weld: Instrument Tee	74
DT-29	Hygienic Clamp Joint: Instrument Tee	74
DT-30	Automatic Tube Weld: Cap	74
DT-V-1	Hygienic Clamp Joint: Weir Style Diaphragm Valve	75
MJ-1	Acceptance Criteria for Welds on Pressure Vessels and Tanks	79
MJ-2	Acceptance Criteria for Welds on Pipe	80
MJ-3	Acceptance Criteria for Welds on Tube	81
MJ-4	Acceptance Criteria for Tube-Attachment Welds	83
SF-1	Acceptance Criteria for Stainless Steel and Higher Alloy Mechanically Polished Product Contact Surface Finishes	87
SF-2	Acceptance Criteria for Mechanically Polished and Electropolished Product Surface Finishes	88
SF-3	R_a Readings for Product Contact Surfaces	88
SG-1	Common Rotary Seal Materials for Biochemical and Sterile Service	96
PM-1	Size Comparison of Common Thermoplastic Sizing Standards	105
Nonmandatory Appendices		
A	Commentary: Slag	111
B	Material Examination Log and Weld Log	112
C	Slope Measurement	117
Index		118